




GABRE Antibody

Product Code	CSB-PA004692
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78334
Immunogen	Synthesized peptide derived from the Internal region of Human GABAA Rε.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GABRE; Gamma-aminobutyric acid receptor subunit epsilon; GABA(A) receptor subunit epsilon
Immunogen Species	Homo sapiens (Human)
Target Names	GABRE

Image



Western Blot analysis of Jurkat cells using GABAA Rε Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.