




NUP88 Antibody

Product Code	CSB-PA004698
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99567
Immunogen	Synthesized peptide derived from the Internal region of Human Nup88.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NUP88; Nuclear pore complex protein Nup88; 88 kDa nucleoporin; Nucleoporin Nup88
Immunogen Species	Homo sapiens (Human)
Target Names	NUP88

Image



Western Blot analysis of K562 cells using Nup88 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.