




MPHOSPH10 Antibody

Product Code	CSB-PA004713
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00566
Immunogen	Synthesized peptide derived from the Internal region of Human MPP10.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MPHOSPH10; MPP10; U3 small nucleolar ribonucleoprotein protein MPP10; M phase phosphoprotein 10
Immunogen Species	Homo sapiens (Human)
Target Names	MPHOSPH10

Image



Western Blot analysis of Jurkat cells using MPHOSPH10 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.