




UBE2A/UBE2B Antibody

Product Code	CSB-PA004716
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49459/P63146
Immunogen	Synthesized peptide derived from the C-terminal region of Human UBE2A/B.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	UBE2A; RAD6A; Ubiquitin-conjugating enzyme E2 A; RAD6 homolog A; HR6A; hHR6A; Ubiquitin carrier protein A; Ubiquitin-protein ligase A; UBE2B; RAD6B; Ubiquitin-conjugating enzyme E2 B;
Immunogen Species	Homo sapiens (Human)
Target Names	UBE2A/UBE2B

Image



Western Blot analysis of Jurkat cells using UBE2A/B Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.