




TRIM68 Antibody

Product Code	CSB-PA004718
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6AZZ1
Immunogen	Synthesized peptide derived from the Internal region of Human SS-56.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TRIM68; GC109; RNF137; SS56; E3 ubiquitin-protein ligase TRIM68; RING finger protein 137; SSA protein SS-56; SS-56; Tripartite motif-containing protein 68
Immunogen Species	Homo sapiens (Human)
Target Names	TRIM68

Image



Western Blot analysis of Jurkat K562 cells using SS-56 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.