




SND1 Antibody

Product Code	CSB-PA004724
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7KZF4
Immunogen	Synthesized peptide derived from the Internal region of Human TudorSN.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SND1; TDRD11; Staphylococcal nuclease domain-containing protein 1; 100 kDa coactivator; EBNA2 coactivator p100; Tudor domain-containing protein 11; p100 co-activator
Immunogen Species	Homo sapiens (Human)
Target Names	SND1

Image



Western Blot analysis of Jurkat cells using TudorSN Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.