




DPP7 Antibody

Product Code	CSB-PA004836
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHL4
Immunogen	Synthesized peptide derived from the C-terminal region of Human DPP7.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DPP7; DPP2; QPP; Dipeptidyl peptidase 2; Dipeptidyl aminopeptidase II; Dipeptidyl peptidase 7; Dipeptidyl peptidase II; DPP II; Quiescent cell proline dipeptidase
Immunogen Species	Homo sapiens (Human)
Target Names	DPP7

Image



Western Blot analysis of Jurkat cells using DPP7 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.