




CA4 Antibody

Product Code	CSB-PA004851
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22748
Immunogen	Synthesized peptide derived from the Internal region of Human CA IV.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CA4; Carbonic anhydrase 4; Carbonate dehydratase IV; Carbonic anhydrase IV; CA-IV
Immunogen Species	Homo sapiens (Human)
Target Names	CA4

Image



Western Blot analysis of rat muscle cells using CA IV Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.