




TXNDC17 Antibody

Product Code	CSB-PA004860
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRA2
Immunogen	Synthesized peptide derived from the Internal region of Human TRP14.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TXNDC17; TXNL5; Thioredoxin domain-containing protein 17; 14 kDa thioredoxin-related protein; TRP14; Protein 42-9-9; Thioredoxin-like protein 5
Immunogen Species	Homo sapiens (Human)
Target Names	TXNDC17

Image



Western Blot analysis of HL-60 Raji K562 SW480 A549 22RV-1 cells using TRP14 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.