




PREP Antibody

Product Code	CSB-PA004861
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6D4
Immunogen	Synthesized peptide derived from the N-terminal region of Human PREP.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PREP; PEP; Prolyl endopeptidase; PE; Post-proline cleaving enzyme
Immunogen Species	Homo sapiens (Human)
Target Names	PREP

Image



Western Blot analysis of Jurkat cells using PREP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.