




NR1H3 Antibody

Product Code	CSB-PA004906
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13133
Immunogen	Synthesized peptide derived from the Internal region of Human LXR α .
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NR1H3; LXRA; Oxysterols receptor LXR-alpha; Liver X receptor alpha; Nuclear receptor subfamily 1 group H member 3
Immunogen Species	Homo sapiens (Human)
Target Names	NR1H3

Image



Western Blot analysis of HeLa cells using LXR α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.