




PPID Antibody

Product Code	CSB-PA004907
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08752
Immunogen	Synthesized peptide derived from the C-terminal region of Human Cyclophilin D.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPID; CYP40; CYPD; Peptidyl-prolyl cis-trans isomerase D; PPlase D; 40 kDa peptidyl-prolyl cis-trans isomerase; Cyclophilin-40; CYP-40; Cyclophilin-related protein; Rotamase D
Immunogen Species	Homo sapiens (Human)
Target Names	PPID

Image



Western Blot analysis of HeLa cells using Cyclophilin D Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.