




NCAM1 Antibody

Product Code	CSB-PA004909
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13591
Immunogen	Synthesized peptide derived from the Internal region of Human CD56.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NCAM1; NCAM; Neural cell adhesion molecule 1; N-CAM-1; NCAM-1; CD56
Immunogen Species	Homo sapiens (Human)
Target Names	NCAM1

Image



Western Blot analysis of H460 MCF7 K562 BT474 mouse brain cells using CD56 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.