




CFLAR Antibody

Product Code	CSB-PA004981
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15519
Immunogen	Synthesized peptide derived from the N-terminal region of Human c-FLIP.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CFLAR; CASH; CASP8AP1; CLARP; MRIT; CASP8 and FADD-like apoptosis regulator; Caspase homolog; CASH; Caspase-eight-related protein; Casper; Caspase-like apoptosis regulatory protein;
Immunogen Species	Homo sapiens (Human)
Target Names	CFLAR

Image



Western Blot analysis of SW480 SKOV3 HepG2 cells using c-FLIP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.