




AMY1A Antibody

Product Code	CSB-PA005016
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04745
Immunogen	Synthesized peptide derived from the N-terminal region of Human Amylase.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AMY1A; AMY1; AMY1B; AMY1; AMY1C; AMY1; Alpha-amylase 1; 1,4-alpha-D-glucan glucanohydrolase 1; Salivary alpha-amylase; AMY2A; Pancreatic alpha-amylase; PA;
Immunogen Species	Homo sapiens (Human)
Target Names	AMY1A

Image



Western Blot analysis of HeLa cells using Amylase Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.