




NTRK3 Antibody

Product Code	CSB-PA005017
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16288
Immunogen	Synthesized peptide derived from the Internal region of Human Trk C.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NTRK3; TRKC; NT-3 growth factor receptor; GP145-TrkC; Trk-C; Neurotrophic tyrosine kinase receptor type 3; TrkC tyrosine kinase
Immunogen Species	Homo sapiens (Human)
Target Names	NTRK3

Image



Western Blot analysis of HeLa cells using Trk C Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.