




ANGPTL4 Antibody

Product Code	CSB-PA005044
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY76
Immunogen	Synthesized peptide derived from the Internal region of Human Angptl4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ANGPTL4; ARP4; HFARP; PGAR; PP1158; PSEC0166; Angiopoietin-related protein 4; Angiopoietin-like protein 4; Hepatic fibrinogen/angiopoietin-related protein; HFARP
Immunogen Species	Homo sapiens (Human)
Target Names	ANGPTL4

Image



Western Blot analysis of MCF7 cells using Angptl4 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.