




MDK Antibody

Product Code	CSB-PA005049
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21741
Immunogen	Synthesized peptide derived from the Internal region of Human MK.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MDK; MK1; NEGF2; Midkine; MK; Amphiregulin-associated protein; ARAP; Midgestation and kidney protein; Neurite outgrowth-promoting factor 2; Neurite outgrowth-promoting protein
Immunogen Species	Homo sapiens (Human)
Target Names	MDK

Image



Western Blot analysis of SW480 MCF7 PC-3 Raji K562 mouse eye mouse kidney cells using MK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.