




CDC34 Antibody

Product Code	CSB-PA005107
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49427
Immunogen	Synthesized peptide derived from the Internal region of Human Cdc34.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDC34; UBCH3; UBE2R1; Ubiquitin-conjugating enzyme E2 R1; Ubiquitin-conjugating enzyme E2-32 kDa complementing; Ubiquitin-conjugating enzyme E2-CDC34; Ubiquitin-protein ligase R1
Immunogen Species	Homo sapiens (Human)
Target Names	CDC34

Image



Western Blot analysis of K562 PC-3 mouse liver cells using Cdc34 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.