




FGF5 Antibody

Product Code	CSB-PA005108
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12034
Immunogen	Synthesized peptide derived from the C-terminal region of Human FGF-5.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FGF5; Fibroblast growth factor 5; FGF-5; Heparin-binding growth factor 5; HBGF-5; Smag-82
Immunogen Species	Homo sapiens (Human)
Target Names	FGF5

Image



Western Blot analysis of mouse liver cells using FGF-5 Polyclonal Antibody