




PPIA Antibody

Product Code	CSB-PA005117
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62937
Immunogen	Synthesized peptide derived from the Internal region of Human CyPA.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPIA; CYPA; Peptidyl-prolyl cis-trans isomerase A; PPIase A; Cyclophilin A; Cyclosporin A-binding protein; Rotamase A
Immunogen Species	Homo sapiens (Human)
Target Names	PPIA

Image



Western Blot analysis of 293 HepG2 A549 SHSY5Y PC-12 BXPC-3 Jurkat mouse kidney mouse lung cells using CyPA Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.