




DBI Antibody

Product Code	CSB-PA005119
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07108
Immunogen	Synthesized peptide derived from the C-terminal region of Human ACBP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DBI; Acyl-CoA-binding protein; ACBP; Diazepam-binding inhibitor; DBI; Endozepine; EP
Immunogen Species	Homo sapiens (Human)
Target Names	DBI

Image



Western Blot analysis of MCF7 U251 mouse liver mouse kidney cells using ACBP Polyclonal Antibody