




ANXA2 Antibody

Product Code	CSB-PA005132
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07355
Immunogen	Synthesized peptide derived from the Internal region of Human Annexin II.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ANXA2; ANX2; ANX2L4; CAL1H; LPC2D; Annexin A2; Annexin II; Annexin-2; Calpactin I heavy chain; Calpactin-1 heavy chain; Chromobindin-8; Lipocortin II; Placental anticoagulant protein IV; PAP-IV; Protein I; p36
Immunogen Species	Homo sapiens (Human)
Target Names	ANXA2

Image



Western Blot analysis of HeLa HepG2 cells using Annexin II Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.