




FSCN1 Antibody

Product Code	CSB-PA005142
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16658
Immunogen	Synthesized peptide derived from the Internal region of Human Fascin 1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FSCN1; FAN1; HSN; SNL; Fascin; 55 kDa actin-bundling protein; Singed-like protein; p55
Immunogen Species	Homo sapiens (Human)
Target Names	FSCN1

Image



Western Blot analysis of A549 HuvEc HT29 K562 L929 cells using Fascin 1 Polyclonal Antibody