




IL5RA Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA005144 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q01344 |
| Immunogen | Synthesized peptide derived from the N-terminal region of Human IL-5R α . |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | IL5RA; IL5R; Interleukin-5 receptor subunit alpha; IL-5 receptor subunit alpha; IL-5R subunit alpha; IL-5R-alpha; IL-5RA; CDw125; CD125 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IL5RA |

Image



Western Blot analysis of K562 cells using IL-5R α Polyclonal Antibody