




IL5RA Antibody

Product Code	CSB-PA005144
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01344
Immunogen	Synthesized peptide derived from the N-terminal region of Human IL-5R α .
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL5RA; IL5R; Interleukin-5 receptor subunit alpha; IL-5 receptor subunit alpha; IL-5R subunit alpha; IL-5R-alpha; IL-5RA; CDw125; CD125
Immunogen Species	Homo sapiens (Human)
Target Names	IL5RA

Image



Western Blot analysis of K562 cells using IL-5R α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.