




PDIA3 Antibody

Product Code	CSB-PA005150
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30101
Immunogen	Synthesized peptide derived from the Internal region of Human ERp57.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PDIA3; ERP57; ERP60; GRP58; Protein disulfide-isomerase A3; 58 kDa glucose-regulated protein; 58 kDa microsomal protein; p58; Disulfide isomerase ER-60;
Immunogen Species	Homo sapiens (Human)
Target Names	PDIA3

Image



Western Blot analysis of A549 HT29 K562 cells using ERp57 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.