




DUSP3 Antibody

Product Code	CSB-PA005173
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51452
Immunogen	Synthesized peptide derived from the C-terminal region of Human VHR.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DUSP3; VHR; Dual specificity protein phosphatase 3; Dual specificity protein phosphatase VHR; Vaccinia H1-related phosphatase; VHR
Immunogen Species	Homo sapiens (Human)
Target Names	DUSP3

Image



Western Blot analysis of HepG2 cells using VHR Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.