




SNAP23 Antibody

Product Code	CSB-PA005176
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00161
Immunogen	Synthesized peptide derived from the C-terminal region of Human SNAP 23.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SNAP23; Synaptosomal-associated protein 23; SNAP-23; Vesicle-membrane fusion protein SNAP-23
Immunogen Species	Homo sapiens (Human)
Target Names	SNAP23

Image



Western Blot analysis of HepG2 cells using SNAP 23 Polyclonal Antibody