




SNAP23 Antibody

Product Code	CSB-PA005176
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00161
Immunogen	Synthesized peptide derived from the C-terminal region of Human SNAP 23.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SNAP23; Synaptosomal-associated protein 23; SNAP-23; Vesicle-membrane fusion protein SNAP-23
Immunogen Species	Homo sapiens (Human)
Target Names	SNAP23

Image



Western Blot analysis of HepG2 cells using SNAP 23 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.