




HEXA Antibody

Product Code	CSB-PA005178
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06865
Immunogen	Synthesized peptide derived from the Internal region of Human HEXA.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HEXA; Beta-hexosaminidase subunit alpha; Beta-N-acetylhexosaminidase subunit alpha; Hexosaminidase subunit A; N-acetyl-beta-glucosaminidase subunit alpha
Immunogen Species	Homo sapiens (Human)
Target Names	HEXA

Image



Western Blot analysis of HepG2 cells using HEXA Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.