





TF Antibody

Product Code	CSB-PA005189
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02787
Immunogen	Synthesized peptide derived from the C-terminal region of Human Transferrin.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TF; Serotransferrin; Transferrin; Beta-1 metal-binding globulin; Siderophilin
Immunogen Species	Homo sapiens (Human)
Target Names	TF

Image



Western Blot analysis of K562 cells using Transferrin Polyclonal Antibody



Western Blot analysis of HELA using TF Polyclonal Antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.