




FCGR2A Antibody

Product Code	CSB-PA005203
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12318
Immunogen	Synthesized peptide derived from the Internal region of Human CD32-A.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FCGR2A; CD32; FCG2; FCGR2A1; IGFR2; Low affinity immunoglobulin gamma Fc region receptor II-a; IgG Fc receptor II-a; CDw32; Fc-gamma RII-a; Fc-gamma-RIIa; FcRII-a; CD32
Immunogen Species	Homo sapiens (Human)
Target Names	FCGR2A

Image



Western Blot analysis of K562 cells using CD32-A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.