




CD59 Antibody

Product Code	CSB-PA005228
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13987
Immunogen	Synthesized peptide derived from the Internal region of Human CD59.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD59; MIC11; MIN1; MIN2; MIN3; MSK21; CD59 glycoprotein; 1F5 antigen; 20 kDa homologous restriction factor; HRF-20; HRF20; MAC-inhibitory protein; MAC-IP;
Immunogen Species	Homo sapiens (Human)
Target Names	CD59

Image



Western Blot analysis of K562 cells using CD59 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.