




CD79B Antibody

Product Code	CSB-PA005233
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40259
Immunogen	Synthesized peptide derived from the Internal region of Human CD79b.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD79B; B29; IGB; B-cell antigen receptor complex-associated protein beta chain; B-cell-specific glycoprotein B29; Ig-beta; Immunoglobulin-associated B29 protein; CD79b
Immunogen Species	Homo sapiens (Human)
Target Names	CD79B

Image



Western Blot analysis of MCF7 L929 K562 cells using CD79b Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.