




CD84 Antibody

Product Code	CSB-PA005235
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIB8
Immunogen	Synthesized peptide derived from the Internal region of Human CD84.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD84; SLAMF5; SLAM family member 5; Cell surface antigen MAX.3; Hly9-beta; Leukocyte differentiation antigen CD84; Signaling lymphocytic activation molecule 5; CD84
Immunogen Species	Homo sapiens (Human)
Target Names	CD84

Image



Western Blot analysis of K562 L929 cells using CD84 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.