




SEMA4D Antibody

Product Code	CSB-PA005238
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92854
Immunogen	Synthesized peptide derived from the N-terminal region of Human CD100.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SEMA4D; C9orf164; CD100; SEMAJ; Semaphorin-4D; A8; BB18; GR3; CD100
Immunogen Species	Homo sapiens (Human)
Target Names	SEMA4D

Image



Western Blot analysis of MCF7 L929 cells using CD100 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.