




FCRL5 Antibody

Product Code	CSB-PA005263
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RD9
Immunogen	Synthesized peptide derived from the Internal region of Human CD307.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FCRL5; FCRH5; IRTA2; Fc receptor-like protein 5; FcR-like protein 5; FcRL5; BXMAS1; Fc receptor homolog 5; FcRH5; Immune receptor translocation-associated protein 2; CD307e
Immunogen Species	Homo sapiens (Human)
Target Names	FCRL5

Image



Western Blot analysis of HepG2 cells using CD307 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.