





PROCR Antibody

Product Code	CSB-PA005324
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNN8
Immunogen	Synthesized peptide derived from the Internal region of Human CD201.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PROCR; EPCR; Endothelial protein C receptor; Activated protein C receptor; APC receptor; Endothelial cell protein C receptor; CD201
Immunogen Species	Homo sapiens (Human)
Target Names	PROCR

Image



Western Blot analysis of K562 cells using CD201 Polyclonal Antibody



Western blot analysis of Rat-brain lysis using CD201 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.