




OGN Antibody

Product Code	CSB-PA005349
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20774
Immunogen	Synthesized peptide derived from the C-terminal region of Human Osteoglycin.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OGN; OIF; SLRR3A; Mimecan; Osteoglycin; Osteoinductive factor; OIF
Immunogen Species	Homo sapiens (Human)
Target Names	OGN

Image



Western Blot analysis of L929 cells using Osteoglycin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.