




PDGFD Antibody

Product Code	CSB-PA005375
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZP0
Immunogen	Synthesized peptide derived from the C-terminal region of Human PDGF-D.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PDGFD; IEGF; SCDGFB; MSTP036; Platelet-derived growth factor D; PDGF-D; Iris-expressed growth factor; Spinal cord-derived growth factor B; SCDGF-B
Immunogen Species	Homo sapiens (Human)
Target Names	PDGFD

Image



Western Blot analysis of L929 cells using PDGF-D Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.