




BMP8A Antibody

Product Code	CSB-PA005425
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5Y6
Immunogen	Synthesized peptide derived from the Internal region of Human BMP-8A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BMP8A; Bone morphogenetic protein 8A; BMP-8A
Immunogen Species	Homo sapiens (Human)
Target Names	BMP8A

Image



Western Blot analysis of Jurkat cells using BMP-8A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.