




ZKSCAN4 Antibody

Product Code	CSB-PA005456
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969J2
Immunogen	Synthesized peptide derived from the Internal region of Human ZKSCAN4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZKSCAN4; ZNF307; ZNF427; Zinc finger protein with KRAB and SCAN domains 4; P373c6.1; Zinc finger protein 307; Zinc finger protein 427
Immunogen Species	Homo sapiens (Human)
Target Names	ZKSCAN4

Image



Western Blot analysis of MCF7 HepG2 293 mouse brain mouse lung cells using ZKSCAN4 Polyclonal Antibody