




IDH3A Antibody

Product Code	CSB-PA005463
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50213
Immunogen	Synthesized peptide derived from the Internal region of Human IDH3A.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IDH3A; Isocitrate dehydrogenase [NAD] subunit alpha, mitochondrial; Isocitric dehydrogenase subunit alpha; NAD(+)-specific ICDH subunit alpha
Immunogen Species	Homo sapiens (Human)
Target Names	IDH3A

Image



Western Blot analysis of MCF-7 HT29 mouse brain mouse heart mouse colon mouse kidney mouse skeletal muscle cells using IDH3A Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.