




E2F4 Antibody

Product Code	CSB-PA005469
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16254
Immunogen	Synthesized peptide derived from Human E2F-4 around the non-acetylation site of K96.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	E2F4; Transcription factor E2F4; E2F-4
Immunogen Species	Homo sapiens (Human)
Target Names	E2F4

Image



Western Blot analysis of Jurkat cells using E2F-4 Polyclonal Antibody