




FOXA1 Antibody

Product Code	CSB-PA005476
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55317
Immunogen	Synthesized peptide derived from Human HNF-3 $\alpha/\beta/\gamma$ around the non-acetylation site of K265.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXA1; HNF3A; TCF3A; Hepatocyte nuclear factor 3-alpha; HNF-3-alpha; HNF-3A; Forkhead box protein A1; Transcription factor 3A; TCF-3A; FOXA2; HNF3B; TCF3B; Hepatocyte nuclear factor 3-beta; HNF-3-beta; HNF-3B; Forkhead box protein A2;
Immunogen Species	Homo sapiens (Human)
Target Names	FOXA1

Image



Western Blot analysis of A549 cells using HNF-3 $\alpha/\beta/\gamma$ Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.