




TP73 Antibody

Product Code	CSB-PA005527
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15350
Immunogen	Synthesized peptide derived from Human p73 around the non-acetylation site of K327.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein
Immunogen Species	Homo sapiens (Human)
Target Names	TP73

Image



Western Blot analysis of A549 cells using p73 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.