




IL2RG Antibody

Product Code	CSB-PA005615
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31785
Immunogen	Synthesized peptide derived from the Internal region of Human IL-2R γ .
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IL2RG; Cytokine receptor common subunit gamma; Interleukin-2 receptor subunit gamma; IL-2 receptor subunit gamma; IL-2R subunit gamma; IL-2RG; gammaC; p64; CD132
Immunogen Species	Homo sapiens (Human)
Target Names	IL2RG

Image



Western Blot analysis of K562 cells using IL-2R γ Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.