




NLRP3 Antibody

Product Code	CSB-PA005618
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P20
Immunogen	Synthesized peptide derived from the Internal region of Human Cryopyrin.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NLRP3; C1orf7; CIAS1; NALP3; PYPAF1; NACHT, LRR and PYD domains-containing protein 3; Angiotensin/vasopressin receptor AII/AVP-like; Caterpillar protein 1.1
Immunogen Species	Homo sapiens (Human)
Target Names	NLRP3

Image



Western Blot analysis of K562 cells using Cryopyrin Polyclonal Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.