





CYBB Antibody

Product Code	CSB-PA005619
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04839
Immunogen	Synthesized peptide derived from the Internal region of Human gp91-phox.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CYBB; NOX2; Cytochrome b-245 heavy chain; CGD91-phox; Cytochrome b(558) subunit beta; Cytochrome b558 subunit beta; Heme-binding membrane glycoprotein gp91phox; NADPH oxidase 2
Immunogen Species	Homo sapiens (Human)
Target Names	CYBB

Image



Western Blot analysis of K562 cells using gp91-phox Polyclonal Antibody.



Western Blot analysis of HEPG2 using gp91-phox Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.